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Stand der Information:
02/2012

Instructions for use

Area of application: Crowns and bridges technique according DIN EN ISO 22674, type 4 alloys.

Composition: (in % by mass)	Technical data (guidelines)
Co 61.0	Density g/cm ³ : 8.3
Cr 28.0	Vickers hardness VH 10: 320
W 8.5	WAK (20-600 °C): 14.5 µm/mK
Si 1.65	Melting range °C: 1.390 – 1.415
	Modulus of elasticity Gpa: 190
Other constituents: Mn-C-Fe	Elongation at rupture A ₅ %: 10

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

CoCr BioStar is a precious free, chromium cobalt-based alloy for the dental application used in dental milling machines (CAD-CAM). It does not contain any nickel, beryllium or gallium. One of the remarkable features is the high corrosion resistance and biocompatibility. CoCr BioStar is suitable for soldering. Its low hardness of only 285 VH 10 allows CoCr BioStar to be easily milled.

Detach of the Framework:

Detach the framework and single parts with suitable hard metal mills or cutting wheels, then slur/grind and finally degrease with an adequate agent.

Firing porcelain:

Standard high sintering porcelains can be used. Follow the relevant manufacturer's instructions for use.

- Clean the framework with a steam-jet or with distilled water in an ultrasonic cleaner.
- An oxide firing (atmospheric at 980 °C, 10 minutes) is not absolutely mandatory but could have advantages.
- After the oxide firing blast always the oxide layer with fresh aluminum oxides 110 µm and clean the surface with a steam-jet.
- The use of a bonder is positive for the metal-ceramic bonding. Wash and opaque firing according the manufacturer's instructions for use.
- We recommend a slow cooling depending to the size of the construction according the manufacturer's instructions for use.

Soldering and welding:

Presolder using SILADENT CrCo solder and the appropriate high-fusing flux or a high-fusing solder for precious bonding alloys. Laser weld using a SILADENT CrCo laser welding rod (REF 102806).

Safety hint:

Metal dust is hazardous to health. For finishing and sandblasting use a suitable extraction system and/or face mask e.g. type FFP3-EN 149-2001!

Contraindications, side effects and interactions with other dental alloys:

In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. If lasting, electrochemically induced, local dysaesthesia arises from contact with other alloys, the restorations must be replaced with other materials.

Further questions: Contact SILADENT technical services
(Tel.: +49 (0) 53 21 – 37 79 25 / 26)
or our sales representatives.

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SILADENT Dr. Böhme & Schöps GmbH